

OB/GYN PRODUCTS





HSG Catheter Set

HSG Catheter Sets are useful in evaluating the causes of abnormal bleeding, menstrual disorders, recurrent miscarriage, and unexplained infertility. They are latex-free and available with or without a silicone balloon.



PRODUCT CODE	CATHETER LENGTH	DIMENSION (OD/MM)	TYPE
GYN05S	300	5FR (1.75)	STIFF
GYN07S	300	7FR (2.2)	STIFF
GYN05F	300	5FR (1.75)	FLEXIBLE
GYN07F	300	7FR (2.2)	FLEXIBLE

Bellota Sonohysterography Catheter

The catheter is to be used for sonohysterography procedure that aid in the analysis of uterine pathology. A placement sheath eliminates the need for tenaculum. This catheter is designed for single patient use.



PRODUCT CODE	CATHETER LENGTH	DIMENSION (OD/MM)	SLEEVE LENGTH	TYPE
GYN205WO	350	1.8	205	WITHOUT CERVICAL CONE
GYN125W	296	1.8	125	WITH CERVICAL CONE
GYN155W	296	1.8	155	WITH CERVICAL CONE



- **Rigid Curette**
- The suction Curette is used for uterine evacuation and curettage in voluntary pregnancy termination or following miscarriage.



PRODUCT CODE	CATHETER LENGTH (MM)	DIMENSION (OD/MM)	TYPE
GYN16S	190	6	STRAIGHT
GYN17S	190	7	STRAIGHT
GYN18S	190	8	STRAIGHT
GYN19S	190	9	STRAIGHT
GYN10S	190	10	STRAIGHT
GYN11S	190	11	STRAIGHT
GYN12C	190	12	STRAIGHT
GYN16C	190	6	CURVED
GYN17C	190	7	CURVED
GYN18C	190	8	CURVED
GYN19C	190	9	CURVED
GYN10C	190	10	CURVED
GYN11C	190	11	CURVED
GYN12C	190	12	CURVED

Word/Barthlin Catheter

The Word Bartholin Catheter set contains a silicone inflatable catheter, scalpel, syringe and needle, all needed in the treatment for Bartholin gland abscesses. Uterine



PRODUCT CODE	CATHETER LENGTH (MM)	DIMENSION (OD/MM)
GYN555	55	5



Uterine Injector

The Uterine Injectors are intended to be used during hysterosonography,hydrotubation, hysterosalpingogram, and salpingoplasties.



PRODUCT CODE	CATHETER LENGTH (MM)	DIMENSION (OD/MM)
GYN555	55	5

Uterine Manipulator

The Uterine Manipulators are designed for Laparoscopic Tubal Occlusion, Diagnostic Laparoscopy, Minilaparotomy, Fertility Exams, and Salpingoplastic procedures where manipulation of the uterus is required.



PRODUCT CODE	CATHETER LENGTH (MM)	DIMENSION (OD/MM)
GYN330	330	4.93
GYN350	350	4.93



Stiff Curette and Flexible Curette

The Stiff and Flexible curettes are used to obtain samples of mucosal lining for histological study in cancer screening, endometrial dating, or bacterial culturing.

PRODUCT CODE	CATHETER LENGTH (MM)	DIMENSION (OD/MM)	TYPE
GYN02430	265	3.0	STIFF
GYN02530	265	3.0	FLEXIBE
GYN02636	265	3.6	EASY



Suction Curette with Syringe

The curette with syringe is used to obtain endometrial tissue lining and superficial layers of the uterine endometrial wall for histological study. Comes with a twist and lock syringe and Randall cutting edge.

PRODUCT CODE	CATHETER LENGTH (MM)	DIMENSION (OD/MM)
GYN223	22	3



**LOKENATH PHARMACEUTICALS PVT.LTD
ASSOCIATED WITH
BIOTRONIX HEALTHCARE INDUSTRIES INC.**

Miami, Florida USA

**o/5, Durgapur Colony, New Alipore,
Kolkata-700053, W.B. INDIA**

Ph: 91-33-24889165

Mobile: 91-9903883464/91-9432643431

E-mail: loknathpharma.com

Web: www.loknathgroup.com