

# RESPIRATORY CATALOG


**LOKENATH PHARMACEUTICALS  
ASSOCIATED WITH  
BIOTRONICS HEALTHCARE, U.S.A**



**Lokenath Group  
(INDIA)**

[www.lokenathgroup.com](http://www.lokenathgroup.com)

## Respiratory Catalog



Our associated company Biotronix Healthcare is an Industry leader of respiratory medical supplies, devices and products. We manufacture, supply and directly support our full line of products that includes nebulizers, oxygen masks, breathing circuits, open and closed suction catheters, yankauers, suction canisters, suction tubing and more. Our heat and moisture exchangers (HMEs) provide heat and humidity that result in patient comfort and safety, without condensation.





## Aerosol Mask Products

▪Aerosol Mask



- Designed for patient comfort
- Clear, soft vinyl
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing
- Can connect 6ml and 20ml nebulizer

-- Does not include tubing or nebulizer

ITEM NO.	SIZE	QTY/BOX
RE-334001	Adult standard	50
RE-334003	Pediatric standard	50

▪Aerosol Mask, Elongated



- Under the chin designed for patient comfort
- Clear, soft vinyl
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing
- Can connect 6ml and 20ml nebulizer
- Does not include tubing or nebulizer

ITEM NO.	SIZE	QTY/BOX
RE-334002	Adult elongated	50
RE-334004	Pediatric elongated	50

▪Tracheostomy Mask



- Designed for patient comfort
- Helix connector for 360° swivel
- 22mm outer connector
- Tracheostomy and laryngectomy therapy
- Does not come with tubing or nebulizer

ITEM NO.	SIZE	QTY/BOX
RE-335001	Adult	50
RE-335002	Pediatric	50



## Aerosol Mask Products

- Nebulizers



- 6ml and 20ml nebulizer
- Anti spill jar with screw cap
- Consistent performance and delivery
- Optimal particle size

ITEM NO.	SIZE	QTY/BOX
RE-341006	6ML	50
RE-341020	20ML	50

Nebulizer with mouth piece



- Designed for patient comfort
- Clear, soft vinyl
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing
- Can connect 6ml and 20ml nebulizer
- Includes 7ft oxygen tubing

ITEM NO.	SIZE	QTY/BOX
RE-334001	Adult standard	50
RE-334003	Pediatric standard	50

- Nebulizer Kit with Elongated Aerosol Mask



- Under the chin designed for patient comfort
- Clear, soft vinyl
- Extremely lightweight and flexible
- Adjustable nose clip
- Accepts 22mm ID aerosol tubing
- Can connect 6ml and 20ml nebulizer
- Includes 7ft oxygen tubing

ITEM NO.	SIZE	QTY/BOX
RE-343001	6ml Nebulizer jar (adult standard)	50
RE-343002	6ml Nebulizer jar (adult elongated)	50
RE-343003	6ml Nebulizer jar (pediatric standard)	50
RE-343004	6ml Nebulizer jar (pediatric elongated)	50
RE-343005	20ml Nebulizer jar (adult standard)	50
RE-343006	20ml Nebulizer jar (adult elongated)	50
RE-343007	20ml Nebulizer jar (pediatric standard)	50
RE-343008	20ml Nebulizer jar (pediatric elongated)	50



## Airways Products

### Guedell Airways



Nontoxic polyethylene construction  
Semi rigid, enclosed centre channel  
Color coded  
Individually packaged  
Latex Free

ITEM NO.	SIZE(MM)	QTY/CS
AN311040	40	48
AN311050	50	48
AN311060	60	48
AN311070	70	48
AN311080	80	48
AN311090	90	48
AN311100	100	48
AN311110	110	48
AN311120	120	48

### Berman Airways



Semi-rigid polyethylene construction  
Two channel  
Individually packaged  
Non-Toxic  
Latex Free

ITEM NO.	SIZE	QTY/CS
AN312040	40	48
AN312050	50	48
AN312060	60	48
AN312070	70	48
AN312080	80	48
AN312090	90	48
AN312100	100	48
AN312110	110	48





## Nasopharyngeal Airways




- Made from soft PVC
- Rounded tip for gentle insertion
- Unique design for secure placement
- Flexible for excellent patient comfort
- Available in a full range of sizes including pediatric.
- Latex Free

ITEM NO.	SIZE	QTY/BOX
AN321025	2.5	10
AN321030	3.0	10
AN321035	3.5	10
AN321040	4.0	10
AN321045	4.5	10
AN321050	5.0	10
AN321055	5.5	10
AN321060	6.0	10
AN321065	6.5	10
AN321070	7.0	10
AN321075	7.5	10
AN321080	8.0	10
AN321085	8.5	10
AN321090	9.0	10



## Breathing Circuits


Bio-Flow Circuit With Expandable Limb, Gas Sampling Line & Breathing Bag



Unique customization capability provides an endless number of circuit configurations  
 Better tube compliance  
 Lower kink resistance  
 Naturally warmed inspired gases  
 Affordable


	LENGTH	CODE	QTY/CASE
Adult	180CM	AN.880582.180.4	50
Pediatric	140CM	AN.880581.140.4	50

▪ Bio-Flow Circuit



	LENGTH	CODE	QTY/CASE
Adult	180CM	AN.880582.180.3	50
Pediatric	140CM	AN.880581.140.3	50

▪ Bio-Flow Circuit With Expandable Limb




	LENGTH	CODE	QTY/CASE
Adult	180CM	AN.880582.180.3L	50
Adult	140CM	AN.880581.140.3L	50




## Respiratory Catalog

### Bio-Flow Circuits With Expandable Limb, Bacterial Viral Filter, Breathing Bag




	LENGTH	CODE	QTY/CASE
Adult	180CM	AN.880582.180.4	50
Pediatric	140CM	AN.880581.140.4	50

### Ventilation Circuits With TWO Water Traps, LIMB (SMOOTH BORE)



	LENGTH	CODE	QTY/CASE
Adult	180CM	RE.880515-1	40
Adult	240CM	RE.880515.24.3	40
Pediatric	150CM	RE.880505-1-1	40

### Ventilation Circuits With TWO Water Traps, LIMB (CORRUGATED)




	LENGTH	CODE	QTY/CASE
Adult	180CM	RE.880515-4	40
Adult	240CM	RE.880515.240.C	40
Pediatric	150CM	RE.880505-4-1	40






• Ventilation Circuits With ONE Water Traps, Limb (SMOOTH BORE)




	LENGTH	CODE	QTY/CASE
Adult	180CM	RE.880515-2	40
Adult	240CM	RE.880515.240.2	40
Pediatric	150CM	RE.880505-2-1	40

(COR) Ventilation Circuits With ONE Water Traps, Limb (CORRUGATED)



	LENGTH	CODE	QTY/CASE
Adult	180CM	RE.880515-3	40
Adult	240CM	RE.880515.240.1	40
Pediatric	150CM	RE.880505-3-1	40

Ventilation Circuits With ONE Water Traps, Limb, Catheter Mount, Bacterial Filter (CORRUGATED)



	LENGTH	CODE	QTY/CASE
Adult	180CM	RE.880515-3A	40
Adult	240CM	RE.880515.240.1A	40
Pediatric	150CM	RE.880505-3A	40




## Respiratory Catalog

### Ventilation Circuits (SMOOTH BORE)




	LENGTH	CODE	QTY/CASE
Adult	180CM	RE.880501.180.0.V	40
Adult	240CM	RE.880501.240.0.V	40
Pediatric	150CM	RE.880512.150.0.V	40

### Ventilation Circuits (CORRUGATED)



	LENGTH	CODE	QTY/CASE
Adult	120CM	RE.880504.120.0.V	40
Adult	160CM	RE.880504.160.0.V	
Adult	180CM	RE.880504.180.0.V	
Adult	240CM	RE.880504.240.0.V	40
Pediatric	100CM	RE.880505.100.0.V	
Pediatric	120CM	RE.880505.120.0.V	
Pediatric	150CM	RE.880505.150.0.V	40


### PROXIMAL PRESSURE LINE CIRCUITS Proximal Pressure Line Circuit Corrugated



	LENGTH	CODE	QTY/CASE
Adult	160CM	RE.880515-5	40
Pediatric	150CM	RE.880515-5-1	40



Proximal Pressure Line Circuit Smoothbore



	LENGTH	CODE	QTY/CASE
Adult	180CM	RE.880515-6	40
Pediatric	150CM	RE.880515-6-1	40






## Cannula's and Tubing Products



Nasal Cannulas

- Designed for patient comfort - Nontoxic medical grade PVC
- Extremely lightweight and flexible
- Over the ear design
- Crush proof and kink resistant

ITEM NO.	SIZE	QTY/BOX
RE-351001	Adult, 7' Tubing	50
RE-351002	Children, 7' Tubing	50
RE-351003	Infant, 7' Tubing	50
RE-351001.14	Adult, 14' Tubing	50
RE-351001.25	Adult , 25' Tubing	50



Oxygen Catheter

- Constructed of nontoxic PVC
- Suitable for adult, children and infant
- Flexible, smooth and odorless
- Apertures are spaced to provide optimal patient comfort

ITEM NO.	SIZE	QTY/BOX
RE-352008	8	50
RE-352010	10	50
RE-352012	12	50
RE-352014	14	50





Oxygen Catheter

- Designed for patient comfort
- Nontoxic medical grade PVC

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-361007S	Standard	7	50
RE-361014S	Standard	14	50
RE-361025S	Standard	25	50
RE-361007U	Universal	7	50
RE-361014U	Universal	14	50
RE-361025U	Universal	25	50



## Oxygen Masks



Medium-Concentration Masks

- Clear, soft vinyl for patient comfort and visual assessment
- Adjustable nose clip assures comfortable fit
- Complete with 7-ft oxygen supply tubing

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-331001	Standard	Adult, 7' Tubing	50
RE-331005	Standard	Children, 7' Tubing	50
RE-331005I	Standard	Infant, 7' Tubing	50
RE-331001U	Universal	Adult, 7' Tubing	50
RE-331005U	Standard	Children, 7' Tubing	50
RE-331005IU	Standard	Infant, 7' Tubing	50
RE-331001.14U	Standard	Adult, 14' Tubing	50
RE-331001.25	Standard	Adult , 25' Tubing	50



Medium-Concentration Masks, Elongated

- Under-the-chin design for excellent fit on wide range of face sizes
- Clear, soft vinyl for patient comfort and visual assessment
- Adjustable nose clip assures comfortable fit
- Complete with 7-ft oxygen supply tubing

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-331003	Standard	Adult, 7' Tubing	50
RE-331007	Standard	Children, 7' Tubing	50
RE-331007I	Standard	Infant, 7' Tubing	50
RE-331003U	Universal	Adult, 7' Tubing	50
RE-331007U	Universal	Children, 7' Tubing	50
RE-331007IU	Universal	Infant, 7' Tubing	50
RE-331003.14	Standard	Adult, 14' Tubing	50
RE-331003.25	Standard	Adult , 25' Tubing	50





## Respiratory Catalog



- Clear, soft vinyl for patient comfort and visual assessment
- Adaptor swivels to accommodate patient position
- Adjustable nose clip assures comfortable fit
- 750-mL reservoir bag
- Complete with 7-ft oxygen supply tubing

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-361007S	Standard	7	50
RE-361014S	Standard	14	50
RE-361025S	Standard	25	50
RE-361007U	Universal	7	50
RE-361014U	Universal	14	50
RE-361025U	Universal	25	50



- Under-the-chin design for excellent fit on wide range of face sizes
- Clear, soft vinyl for patient comfort and visual assessment
- Adaptor swivels to accommodate patient position
- Adjustable nose clip assures comfortable fit
- 750-mL reservoir bag
- Complete with 7-ft oxygen supply tubing

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-361007S	Standard	7	50
RE-361014S	Standard	14	50
RE-361025S	Standard	25	50
RE-361007U	Universal	7	50
RE-361014U	Universal	14	50
RE-361025U	Universal	25	50



## Respiratory Catalog



- Low-resistance check valve does not contain natural rubber latex, prevents rebreathing and allows exhaled gas to escape
- Adjustable nose clip assures comfortable fit
- 750-mL reservoir bag
- Complete with 7-ft oxygen supply tubing

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-361007S	Standard	7	50
RE-361014S	Standard	14	50
RE-361025S	Standard	25	50
RE-361007U	Universal	7	50
RE-361014U	Universal	14	50
RE-361025U	Universal	25	50



- 3-in-1 mask can be used as a medium-concentration, high-concentration or nonbreathing mask
- Includes mask with valve, nonbreathing bag assembly, 7-ft oxygen supply tubing and oxygen connector
- Adjustable nose clip assures comfortable fit
- 750-mL reservoir bag

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-333001S	Standard	Adult Standard	50
RE-333002S	Standard	Adult Elongated	50
RE-333003S	Standard	Pediatric Standard	50
RE-333004S	Standard	Pediatric Elongated	50
RE-333001U	Universal	Adult Standard	50
RE-332002U	Universal	Adult Elongated	50
RE-333003U	Universal	Pediatric Standard	50
RE-333004U	Universal	Pediatric Elongated	50




## Respiratory Catalog



- Safe, simple delivery of variable oxygen concentrations
- Each mask includes color-coded diluters: green for low concentration, white for medium concentration
- Locking ring secures flow setting
- Includes adaptor for high humidity entrainment
- Complete with 7-ft oxygen tubing

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-332001S	Standard	Adult Standard	50
RE-332002S	Standard	Adult Elongated	50
RE-332003S	Standard	Pediatric Standard	50
RE-332004S	Standard	Pediatric Elongated	50
RE-332001U	Universal	Adult Standard	50
RE-332002U	Universal	Adult Elongated	50
RE-332003U	Universal	Pediatric Standard	50
RE-332004U	Universal	Pediatric Elongated	50



- Clear, soft vinyl for patient comfort and visual assessment
- Helix Connector
- Secure Strap

ITEM NO.	CONNECTOR	SIZE	QTY/BOX
RE-335001S	Standard	Adult	50
RE-335002S	Standard	Pediatric	50
RE-335001U	Universal	Adult	50
RE-335002U	Universal	Pediatric	50

